

Anti-ADCY2 antibody (Multiple Epitope) {Biotin} (STJ507775)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/IP/WB |
| Host / Source | Mouse |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene Symbol | Adcy2 |
| UniProt ID | ADCY2_MOUSE |
| Immunogen Region | Multiple Epitope |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|