

## Adeno-Associated Virus Capsid Protein Blocking Peptide for STJ507794 peptide (STJ507797)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
Reactivity	Virus

### PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	C-Term
Specificity	This blocking peptide is recommended for use in combination with Adeno-Associated Virus Capsid Protein antibody, STJ507794

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---