

Factor VIII Ig Gamma Chain Positive Control for STJ507798 peptide (STJ507802)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	C-Term
Specificity	This is positive control is recommended for use in combination with Factor VIII Ig Gamma Chain antibody STJ507798.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---