

Dystrophin Blocking Peptide for STJ507807 peptide (STJ507810)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1756
Gene Symbol	DMD
UniProt ID	DMD_HUMAN
Immunogen Region	200-250
Specificity	This blocking peptide is recommended for use in combination with Dystrophin antibody, STJ507807
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.