

Dystrophin Positive Control for STJ507808 peptide (STJ507811)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	WB
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1756
Gene Symbol	DMD
UniProt ID	DMD_HUMAN
Specificity	This is positive control is recommended for use in combination with Dystrophin antibody STJ507808.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.