

SATB2 Blocking Peptide for STJ507820 peptide (STJ507823)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	23314
Gene Symbol	SATB2
UniProt ID	SATB2_HUMAN
Immunogen Region	400-450
Specificity	This blocking peptide is recommended for use in combination with SATB2 antibody, STJ507820
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.