

Dopamine beta hydroxylase Positive Control for STJ507829 peptide (STJ507834)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13166
Gene Symbol	Dbh
UniProt ID	DOPO_MOUSE
Immunogen Region	125-178
Specificity	This is positive control is recommended for use in combination with Dopamine beta hydroxylase antibody STJ507829.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---