

## PGLUT3 Blocking Peptide for STJ507830 peptide (Phospho) (STJ507845)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	450-523
Specificity	This blocking peptide is recommended for use in combination with PGLUT3 antibody, STJ507830

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---