

PLAP Blocking Peptide for STJ507847 peptide (STJ507850)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	250
Gene Symbol	ALPP
UniProt ID	PPB1_HUMAN
Immunogen Region	N-Term
Specificity	This blocking peptide is recommended for use in combination with PLAP antibody, STJ507847

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---