

Anti-DAPK2 antibody (C-Term) (STJ70101)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Pig

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 40-45kDa bands seen in mouse brain extracts. Recommended for use at 1-3µg/ml. Primary incubation was for 1 hour. IHC: Paraffin embedded Human Breast and Thyroid. Recommended
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23604
Gene Symbol	DAPK2
UniProt ID	DAPK2_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	KALHPRRRSSTS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70101 staining (2µg/ml) of mouse brain extracts (RIPA buffer, 35µg total protein per lane). Detected by chemiluminescence.