

Anti-SCAP2/PRAP antibody (C-Term) (STJ70135)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IF/FC
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES



Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:3000. WB: Approx. 26+37kDa bands observed in lysates of cell line K562, and approx. 28+38kDa in preliminary testing of U937 cell lysate (calculated MW of 21.6kDa according to NP_001290397.1 and 4
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8935
Gene Symbol	SKAP2
UniProt ID	SKAP2_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	GLVPKAYIMEMYDI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



HEK293 overexpressing SKAP2 (RC206119) and STJ70135 (0.5 µg/ml) staining of Jurkat lysate (35 µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.