

Anti-UEV2V2/MMS2 antibody (N-Term) (STJ70164)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 20kDa band seen in human brain lysate. Recommended for use at 0.2-1µg/ml.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7336
Gene Symbol	UBE2V2
UniProt ID	UB2V2_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	AVSTGVKVPKRF
Specificity	This antibody is also expected to recognize the readthrough transcript product KUA-UEV (NP_954673.1, GeneID 387522). Also cross-reactivity is expected with UEV2V1 (GeneID 7335).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70164 staining (0.5µg/ml) of human brain lysate (RIPA buffer, 30µg total protein per lane). Primary incubated for 12 hour. Detected by western blot using chemiluminescence.