

## Anti-LIS1/PAFAH1B1 antibody (C-Term) (STJ70179)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | Pep-ELISA/WB/IHC   |
| Host / Source | Goat               |
| Reactivity    | Human/Rat          |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                                                                                                                |
| Concentration       | 0.5 mg/mL                                                                                                                                                                                                                                                 |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                              |
| Purification        | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.                                                                                                                     |
| Dilution Range      | Peptide ELISA: antibody detection limit dilution 1:32000.<br>WB: Approx 40kDa band observed in Human Lung and Brain (Cerebellum) and in Rat Ovary lysates (calculated MW of 46.6kDa according to Human NP_000421.1 and to Rat NP_113951.1). Recommended c |
| Formulation         | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA                                                                                                                                                                    |
| Isotype             | IgG                                                                                                                                                                                                                                                       |
| Storage Instruction | Store at -20°C on receipt and minimise freeze-thaw cycles.                                                                                                                                                                                                |

### TARGET INFORMATION

|                    |                            |
|--------------------|----------------------------|
| Gene ID            | <a href="#">5048</a>       |
| Gene Symbol        | <a href="#">PAFAH1B1</a>   |
| UniProt ID         | <a href="#">LIS1_HUMAN</a> |
| Immunogen Region   | C-Term                     |
| Immunogen Sequence | TGSVDQTVKVWECR             |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70179 (0.3 µg/ml) staining of Rat Ovary lysate (35 µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.