

Anti-RNF31/ZIBRA antibody (C-Term) (STJ70235)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 125kDa band observed in lysates of cell line JurkatT (calculated MW of 120kDa according to AAH12077.1). Recommended concentration: 0.5-1.5µg/ml.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55072
Gene Symbol	RNF31
UniProt ID	RNF31_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	EEVPLGQSIPRRRK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70235 (0.5µg/ml) staining of JURKAT lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.