

## Anti-GOLPH2 antibody (C-Term) (STJ70257)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. IHC: In paraffin embedded Human Prostate shows luminal staining of secretory cells in the glands. Recommended concentration: 3-6µg/ml. Paraffin embedded Human Uterus and Small Intestine. R
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">51280</a>
Gene Symbol	<a href="#">GOLM1</a>
UniProt ID	<a href="#">GOLM1_HUMAN</a>
Immunogen Region	C-Term
Immunogen Sequence	NLLDQREKRNHTL
Specificity	Both variants (NP_057632.2; NP_808800.1) encode the same isoform

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---