

Anti-UBC3B antibody (N-Term) (STJ70292)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IEM/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: No signal obtained yet but low background observed in Human Testis and U937 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field-have any results been repo
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54926
Gene Symbol	UBE2R2
UniProt ID	UB2R2_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	AQQQMTSSQKALM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70292 (2µg/ml) staining of paraffin embedded Human Testis. Steamed antigen retrieval with citrate buffer pH 6, HRP-staining.