

Anti-EVC2/Limbin antibody (C-Term) (STJ70349)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: A customer reported that a consistent band was observed in 293 and human heart at approx 75kDa with 1 hour incubation of primary at 0.5ug/ml. IHC: In paraffin embedded Human Placen
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	132884
Gene Symbol	EVC2
UniProt ID	LBN_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	LNAKKAMRALGMD
Specificity	This antibody is expected to recognise isoform 1 (NP_667338.3) and isoform 2 (NP_001159608.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



STJ70349 (2µg/ml) staining of paraffin embedded Human Placenta. Steamed antigen retrieval with citrate buffer pH 6, HRP-staining.