

Anti-RANBPM/RANBP9 antibody (C-Term) (STJ70395)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC/ICC
Host / Source	Goat
Reactivity	Human/Mouse/Dog/Pig

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: Preliminary experiments gave no signal but low background in whole cell Hela lysates at up to 1 µg/ml. IHC: In paraffin embedded Human Tonsil stains nuclei in epithelial cells of to
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10048
Gene Symbol	RANBP9
UniProt ID	RANB9_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	GSCAFATVEDYLH
Specificity	This product is also expected to recognize the centrosome-associated truncated isoform RANBPM (BAA23216.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70395 (3 µg/ml) staining of paraffin embedded Human Tonsil. Microwaved antigen retrieval with Tris/EDTA buffer pH9, HRP-staining.