

Anti-BAF57/SMARCE1 antibody (C-Term) (STJ70407)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx 50kDa band observed in lysates of cell lines HEPG2, HeLa and Jurkat (calculated MW of 46.6kDa according to NP_003070.3). A 48kDa band was observed in transfected HEK293 transien
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	6605
Gene Symbol	SMARCE1
UniProt ID	SMCE1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	PPTDPIPEDEKKE
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.