

Anti-NUMBL antibody (N-Term) (STJ70424)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB-Trf
Host / Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: In transfected HEK293 transiently expressing Human Numbl and DYKDDDDK tagged) , a band of approx. 85kDa is observed. No bands are observed in mock-transfected HEK293 and the same band
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9253
Gene Symbol	NUMBL
UniProt ID	NUMBL_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	SRSAAASGGPRRPE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human Numbl with C-terminal MYC tag probed with STJ70424 (1ug/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with STJ70424 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.