

Anti-ARAP3 antibody (C-Term) (STJ70429)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IP
Host / Source	Goat
Reactivity	Human/Mouse/Dog/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: This antibody has been successfully used in WB on scirrhous human gastric cancer cell lines. IP: This antibody has been successfully used in IP on scirrhous human gastric cancer cel
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64411
Gene Symbol	ARAP3
UniProt ID	ARAP3_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	CTSSPPSSQPLT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---