

Anti-CENTD2/ARAP1 antibody (C-Term) (STJ70441)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: This antibody showed a 150KDa in NIH3T3 lysates. IHC: Paraffin embedded Human Lung. Recommended concentration: 6-8µg/ml.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116985
Gene Symbol	ARAP1
UniProt ID	ARAP1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	AFTADPLSLLRNV
Specificity	This antibody is expected to recognise isoform a (NP_056057.2) , isoform c (NP_001035207.1) and isoform d (NP_001128662.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---