

Anti-GNRPX/PLEKHJ1 antibody (C-Term) (STJ70448)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: No signal obtained yet but a preliminary experiment showed low background in Human Liver and Molt-4 extract at up to 1µg/ml. However this used our routine western blotting protocol whi
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	55111
Gene Symbol	PLEKHJ1
UniProt ID	PKHJ1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	SEEARFQLSGLQA
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.