

Anti-NUP62 antibody (C-Term) (STJ70481)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA/WB/FC
Host / Source	Goat
Reactivity	Human/Dog/Cow
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx. 70kDa band observed in nuclear lysates of cell lines Jurkat, K562, A431 and HepG2 (calculated MW of 53.3kDa according to NP_036478.2). This molecular weight is routinely observed
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	23636
Gene Symbol	NUP62
UniProt ID	NUP62_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	RKEQERSFRITFD
Specificity	All four variants represent identical protein.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.