

Anti-Munc13-4/UNC13D antibody (C-Term) (STJ70556)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: This antibody has been successfully used in WB on Human: Cetica V et al, J Allergy Clin Immunol. 2015 May;135 (5) :1310-8. PMID: 25312756. IHC: Paraffin embedded Human Spleen. Reco
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	201294
Gene Symbol	UNC13D
UniProt ID	UN13D_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	KQASQHALRPAP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---