

Anti-DKK1 antibody (C-Term) (STJ70585)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IF/IHC/FC
Host / Source	Goat
Reactivity	Human/Mouse/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 38kDa band observed in lysates of cell line MCF7 and in preliminary testing of HepG2 cell lysate and Human Spleen lysate (calculated MW of 28.7kDa according to NP_036374.1). Thi
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22943
Gene Symbol	DKK1
UniProt ID	DKK1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	CDHHQASNSSLHT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70585 (1 µg/ml) staining of Human Bone Marrow lysate (35 µg protein in RIPA buffer). Detected by chemiluminescence.