

Anti-VDP p115 antibody (C-Term) (STJ70802)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000. WB: Approx 110kDa band observed in Human Kidney and Liver lysates (calculated MW of 108kDa according to NP_003706.1). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8615
Gene Symbol	USO1
UniProt ID	USO1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	DESEDPGKDLDH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70802 (1µg/ml) staining of Human Kidney and Liver lysates (35µg protein in RIPA buffer). Detected by chemiluminescence.