

Anti-EPB41L3/DAL1 antibody (C-Term) (STJ70820)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000. WB: Approx 150kDa band observed in Human Brain (Cerebellum) , Mouse Brain and Rat Brain lysates (calculated MW of 121kDa according to NP_036439.2). Recommended concentration: 0.3-1µg/ml. Pr
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23136
Gene Symbol	EPB41L3
UniProt ID	E41L3_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	HKETEITPEDGED

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70820 (0.3 µg/ml) staining of Human cerebellum lysate (35 µg protein in RIPA buffer). Detected by chemiluminescence.