

Anti-ZNF9/CNBP antibody (C-Term) (STJ70911)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IF/FC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: Approx 20kDa band observed in Human Skeletal Muscle lysates and in lysates of cell lines Daudi, MCF7, MOLT-4 and Neuro-2a (calculated MW of 19.5kDa according to Human NP_003409.1 and Mo
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7555
Gene Symbol	CNBP
UniProt ID	CNBP_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	GESGHLARECTIE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ70911 (2µg/ml) staining of Human Muscle lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.