

Anti-ZNF217 antibody (Internal) (STJ70953)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA/IHC/IF/FC
Host / Source	Goat
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. IHC: In paraffin embedded Human Prostate shows nuclear staining in select secretory cells. Recommended concentration: 3-6µg/ml. IF: Strong expression of the protein seen in the nuclei of
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	7764
Gene Symbol	ZNF217
UniProt ID	ZN217_HUMAN
Immunogen Region	Internal
Immunogen Sequence	KNAADDSDADKVNKN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.