

Anti-FLNB antibody (Internal) (STJ71037)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: Preliminary experiments in Human Brain, Mouse Brain and Rat Brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any f
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2317
Gene Symbol	FLNB
UniProt ID	FLNB_HUMAN
Immunogen Region	Internal
Immunogen Sequence	KITDDSRRC SQVK
Specificity	This antibody is expected to recognise all reported isoforms (NP_001448.2; NP_001157789.1; NP_001157790.1; NP_001157791.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.