

Anti-ACYP1-Isoform a antibody (Internal) (STJ71076)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 15kDa band observed in Human Testis lysates (calculated MW of 11.3kDa according to NP_001098.1). Recommended concentration: 0.01-0.03µg/ml. Primary incubation was 1 hour. Preli
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	97
Gene Symbol	ACYP1
UniProt ID	ACYP1_HUMAN
Immunogen Region	Internal
Immunogen Sequence	PKSHIDKANFNNEK
Specificity	This antibody is expected to recognize isoform a (NP_001098.1) only.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ71076 (0.01µg/ml) staining of Human Testis lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.