

Anti-GPSM2 antibody (Internal) (STJ71084)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC/IF
Host / Source	Goat
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. IHC: Paraffin embedded Human Testis. Recommended concentration: 5µg/ml. IF: Strong expression of the protein seen in the cytoplasm/centrosomes of A431 cells. Recommended concentration:
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29899
Gene Symbol	GPSM2
UniProt ID	GPSM2_HUMAN
Immunogen Region	Internal
Immunogen Sequence	NRLKGKKYKTNSSTK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---