

Anti-STAG3-Mouse antibody (Internal) (STJ71189)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 30kDa band in Human Testis lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature f
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	10734
Gene Symbol	STAG3
UniProt ID	STAG3_HUMAN
Immunogen Region	Internal
Immunogen Sequence	KHYNKFYEDYGD
Specificity	The immunizing peptide was designed based on the Mouse protein sequence with one residue of a difference from the human sequence.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.