

Anti-PHLPP2 antibody (C-Term) (STJ71236)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 140kDa band observed in Human Brain (Amygdala, Hippocampus, Substantia Nigra) lysates (calculated MW of 134kDa according to NP_055835.2). Recommended concentration: 0.1-0.3µg/ml
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23035
Gene Symbol	PHLPP2
UniProt ID	PHLP2_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	PHEEDRTEPPEEFD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ71236 (0.1 µg/ml) staining of Human Brain (Amygdala, Hippocampus, Substantia Nigra) lysate (35 µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.