

Anti-Her6 antibody (N-Term) (STJ71294)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	Pep-ELISA
	Host / Source	Goat
	Reactivity	Zebrafish
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.5 mg/mL
	Conjugation	Unconjugated
	Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
	Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the
	Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
	Isotype	IgG
	Storage Instruction	Store at-20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	N-Term
	Immunogen Sequence	KRKIKQTVTSLPHQ
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.