

Anti-SOX8 antibody (247-258) (STJ71374)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: Preliminary experiments in Human Brain (Cerebellum) , Mouse Brain and Rat Brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	30812
Gene Symbol	SOX8
UniProt ID	SOX8_HUMAN
Immunogen Region	247-258
Immunogen Sequence	EGRRPVDSGRQN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.