

Anti-SEN6/SUSP1 antibody (Internal) (STJ71464)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IF/FC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. IF: Strong expression of the protein seen in the nuclei and cytoplasm of U2OS cells. Recommended concentration: 10 µg/ml. FC: Flow cytometric analysis of U2OS cells. Recommended concentration: 10 µg/ml.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26054
Gene Symbol	SEN6
UniProt ID	SEN6_HUMAN
Immunogen Region	Internal
Immunogen Sequence	KPKYEPNPHYHEN
Specificity	This antibody is expected to recognise both reported isoforms (NP_001093879.1; NP_056386.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ71464 (0.5 µg/ml) staining of HeLa cell nuclear lysate (35 µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.