

## Anti-TRIB1 antibody (C-Term) (STJ71535)

| GENERAL INFORMATION    |                                                                                                                                                                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Primary antibodies                                                                                                                                                                                                                                       |
| Applications           | Pep-ELISA/IHC                                                                                                                                                                                                                                            |
| Host / Source          | Goat                                                                                                                                                                                                                                                     |
| Reactivity             | Human/Mouse/Cow                                                                                                                                                                                                                                          |
| PRODUCT PROPERTIES     |                                                                                                                                                                                                                                                          |
| Clonality              | Polyclonal                                                                                                                                                                                                                                               |
| Concentration          | 0.5 mg/mL                                                                                                                                                                                                                                                |
| Conjugation            | Unconjugated                                                                                                                                                                                                                                             |
| Purification           | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.                                                                                                                    |
| Dilution Range         | Peptide ELISA: antibody detection limit dilution 1:8000.<br>WB: Preliminary experiments gave an approx 35kDa band in Human Lymph Node and Peripheral Blood Mononucleocyte lysates after 0.3µg/ml antibody staining. Please note that currently we cannot |
| Formulation            | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA                                                                                                                                                                   |
| Isotype                | IgG                                                                                                                                                                                                                                                      |
| Storage Instruction    | Store at -20°C on receipt and minimise freeze-thaw cycles.                                                                                                                                                                                               |
| TARGET INFORMATION     |                                                                                                                                                                                                                                                          |
| Gene ID                | <a href="#">10221</a>                                                                                                                                                                                                                                    |
| Gene Symbol            | <a href="#">TRIB1</a>                                                                                                                                                                                                                                    |
| UniProt ID             | <a href="#">TRIB1_HUMAN</a>                                                                                                                                                                                                                              |
| Immunogen Region       | C-Term                                                                                                                                                                                                                                                   |
| Immunogen Sequence     | TSDQIVPEYQEDSD                                                                                                                                                                                                                                           |
| ADDITIONAL INFORMATION |                                                                                                                                                                                                                                                          |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                                                                                                  |