

Anti-CHRFAM7A/CHRNA7-FAM7A antibody (N-Term) (STJ71540)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:1000. WB: Preliminary experiments in Human Brain (Amygdala, Hippocampus, Cerebellum) lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89832
Gene Symbol	CHRFAM7A
UniProt ID	CRFM7_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	QKYCIYQHFQFQ
Specificity	This antibody is expected to recognize reported isoform NP_647536.1 only.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---