

Anti-NMNAT3 antibody (C-Term) (STJ71545)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary testing showed a band at approx 30kDa in Human Placenta and Spleen lysate, and approx. 28kDa in Human Bone Marrow lysate after 0.3-0.5µg/ml antibody staining (calculated)
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	349565
Gene Symbol	NMNAT3
UniProt ID	NMNA3_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	GSTWKGKSTQSTE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---