

Anti-DAGLA antibody (C-Term) (STJ71663)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 125kDa band observed in Human Liver lysates (calculated MW of 115kDa according to NP_006124.1). Recommended concentration: 0.3-1µg/ml. Some minor background is detected and is
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	747
Gene Symbol	DAGLA
UniProt ID	DGLA_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	PAKQDELVISAR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ71663 (0.3µg/ml) staining of Human Liver lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.