

## Anti-RCBTB2 antibody (Internal) (STJ71787)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Cow/Pig

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000. WB: Preliminary experiments gave an approx. 30kDa band in Human Liver, Lung and Tonsil lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in th
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1102</a>
Gene Symbol	<a href="#">RCBTB2</a>
UniProt ID	<a href="#">RCBT2_HUMAN</a>
Immunogen Region	Internal
Immunogen Sequence	CEHFRSSLEDNEDD
Specificity	This antibody is not expected to cross-react with RCBTB1.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---