

Anti-PTMA antibody (C-Term) (STJ71788)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human/Rat/Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:1000. WB: Preliminary experiments gave an approx. 38kDa band in lysates of cell line A431 after 0.2µg/ml antibody staining. Please note that currently we cannot find an explanation in the literat
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5757
Gene Symbol	PTMA
UniProt ID	PTMA_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	DVDTKKQKTDEDD
Specificity	This antibody is expected to recognize both reported isoforms (NP_001092755.1; NP_002814.3).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---