

Anti-ZNF228/ZFP112 antibody (Internal) (STJ71810)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: Preliminary experiments gave an approx. 48kDa band in lysates of cell line HepG2 after 2µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	7771
Gene Symbol	ZNF112
UniProt ID	ZN112_HUMAN
Immunogen Region	Internal
Immunogen Sequence	EREDDIENSHLKSQ
Specificity	This antibody is expected to recognize both reported isoforms (NP_001076804.1; NP_037512.3)
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.