

Anti-UNC5C antibody (Internal) (STJ71863)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IF/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Dog

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments in Human and Mouse Brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). IHC: In paraffin embedded Human Cerebr
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8633
Gene Symbol	UNC5C
UniProt ID	UNC5C_HUMAN
Immunogen Region	Internal
Immunogen Sequence	NCTVSEPTGID

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ71863 (50ug/ml) staining of islets in Mouse Pancreas. Citrate antigen retrieval. Detected by immunofluorescence.