

Anti-Klk1-Mouse antibody (C-Term) (STJ71897)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	Pep-ELISA
	Host / Source	Goat
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.5 mg/mL
	Conjugation	Unconjugated
	Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
	Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: Preliminary experiments gave an approx 55kDa band in Mouse Liver, Heart and Testis lysates after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in
	Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
	Isotype	IgG
	Storage Instruction	Store at-20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	C-Term
	Immunogen Sequence	FNTWIRETMAEND
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.