

Anti-MOZ/KAT6A antibody (Internal) (STJ71939)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: Preliminary experiments in lysates of cell line U937 gave no specific signal but low background (at antibody concentration up to 1µg/ml). IHC: Paraffin embedded Human Adrenal Gland
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	7994
Gene Symbol	KAT6A
UniProt ID	KAT6A_HUMAN
Immunogen Region	Internal
Immunogen Sequence	EKIKDKEETELD
Specificity	Reported variants NP_001092882.1, NP_006757.2, NP_001092883.1 represent identical protein.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.