

Anti-NFATC4/NFAT3 antibody (Internal) (STJ71974)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. WB: Preliminary experiments in Human Heart, Jurkat and MOLT-4 lysates gave no specific signal but low background (at antibody concentration up to 0.5µg/ml). We would appreciate any feedback
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4776
Gene Symbol	NFATC4
UniProt ID	NFAC4_HUMAN
Immunogen Region	Internal
Immunogen Sequence	EATVNRLQSNEVTLT
Specificity	This antibody is expected to recognize both isoforms (NP_001129494.1; NP_004545.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---