

Anti-UPF3A/RENT3A antibody (Internal) (STJ72041)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Cow
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: Preliminary experiments gave an approx 48kDa band in lysates of cell lines HeLa and Kelly after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in t
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	65110
Gene Symbol	UPF3A
UniProt ID	REN3A_HUMAN
Immunogen Region	Internal
Immunogen Sequence	KETDKQKKIAEKE
Specificity	This antibody is expected to recognize both isoforms (NP_075387.1; NP_542418.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.